

**MGHS Test Change Implementation Notice**  
Effective immediately or as otherwise noted**DNA ANTIBODY, DOUBLE STRANDED**  
Marquette General Health System**General Info**

Alpha Code	DNA
MGH LIS Test No	828
Schedule	Monday-Friday
Testing Time	1 day
Testing Lab	Marquette General Hospital
QORR Test Code	DNA

**Specimen Info**

Type	1 SST Tube
Volume	1 ml
Temperature	Refrigerate
Preservative	
Collection Info	Refrigerate specimen after collection.

**Components of Profile**

DNA by Evolis	DNAEV
DNA by IFA	DNA.
DNA Titer	DNAT

\*If EIA screen is negative, DNA. And DNAT are then automatically cancelled by the system.

**Methods**

EIA Screen. All positive results are then tested by Indirect Fluorescence Antibody (IFA) using *Crithidia luciliae* Substrate.

**Clinical Utilities**

Useful for determining antibodies to double-stranded DNA. If Positive results are titered.

**CPT Codes**

86255 - screen 86256 - titer

**Reference Range**

Negative

\*The components and methodology of this assay have been updated. An EIA method is first used to screen all specimens; any positive results are then tested by an IFA method.